Cortisol ELISA

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(Westermann J, Demir A, Herbst V. Determination of Cortisol in Saliva and Serum by a Luminescence-Enhanced Enzyme Immunoassay. Clin Lab 2004;50:11-24)

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- 1 Determination of Cortisol in Saliva and Serum by a Luminescence-Enhanced Enzyme Immunoassay (IBL Jürgen Westermann, Anke Demir, Victor Herbst). *Clin.Lab.* 50 11 - 24 (2004)
- 2 Kirschbaum C. et al. Impact of Gender, Menstrual Cycle Phase and Oral Contraceptives on the Activity of the Hypothalamus-Pituitary-Adrenal Axis. *Psychosomatic Medicine* 61 154 – 162 (1999)
- 3 Raff H., J. L. Raff and J. W. Findling. Late-Night Salivary Cortisol as a Screening Test for Cushing's Syndrome. *J Clin Endocrinol Metab* 83 2681 – 2686 (1998)
- 4 Morineau G. et al. Radioimmunoassay of cortisone in serum, urine and saliva to assess the status of the cortisol-cortisone shuttle. *Clin Chem* 43 1397 – 1407 (1997)
- 5 Tunn S. et al. Simultaneous measurement of cortisol in serum and saliva after different forms of cortisol administration. *Clin Chem* 38 1491 – 1494 (1992)
- 6 Shimada M. et al. Determination of salivary cortisol by ELISA and its application to the assessment of the circadian rhythm in children. *Hum Res* 44 213 – 217 (1995)
- 7 Thomopoulos P. et al. Long distance follow up of a patient with intermittent Cushing's disease by use of salivary cortisol measurements. In: Kirschbaum C. et al. (ed.) *Assessment of hormones and drugs in saliva in biobehavioral research*. 1st Edition. Hofgrebe & Huber. Seattle (1992) pp 89 – 92
- 8 Kahn J. P. et al. Applications of salivary cortisol determinations to psychiatric and stress research: stress responses in students during academic examinations. In: Kirschbaum C. et al. (ed.) *Assessment of hormones and drugs in saliva in biobehavioral research*. 1st Edition. Hofgrebe & Huber. Seattle (1992) pp 111 – 127
- 9 Laudet M. H. et al. Salivary Cortisol Measurement: A Practical Approach to Assess Pituitary-Adrenal Function. *J Clin Endocrinol Metab* 66 343 – 348 (1988)
- 10 Umeda T. et al. Use of Saliva for Monitoring Unbound Free Cortisol Levels in Serum. *Clin Chem Acta* 110 245 – 253 (1981)
- 11 Peters J. R. et al. Salivary Cortisol Assays for Assessing Pituitary-Adrenal Reserve. *Clin Endocrinol* 17 583 – 592 (1982)
- 12 Buchanan et al. Circadian regulation of cortisol after hippocampal damage in humans; *BIOL PSYCHIATRY* 2004; 56:651-656
- 13 Hellhammer et al. Effects of soy lecithin phosphatidic acid and phosphatidylserine complex (PAS) on the endocrine and psychological responses to mental stress; *Stress* 2004 Vol. 00 (0), pp. 1-8
- 14 Gaab et al. Psychological determinations of the cortisol stress response: the role of anticipatory cognitive appraisal; *Psychoneuroendocrinology* 2005 30, 599-610
- 15 Fries et al. Attenuation of the hypothalamic-pituitary-adrenal axis responsivity of the Trier social stress test by the benzodiazepine alprazolam; *Psychoneuroendocrinology* 2006 31, 1278-1288
- 16 Miller et al. clinical depression and regulation of the inflammatory response during acute stress; *Psychosomatic Medicine* 2005 67:679-687
- 17 Kuhlmann et al. Effects of oral cortisol treatment in healthy young women on memory retrieval of negative and neutral words; *Neurobiology of learning and memory* 2005 83 158-162
- 18 Li et al. Life-time socio-economic position and cortisol patterns in mid-life; *Psychoneuroendocrinology* 2007 32, 824-833
- 19 Wolf et al. No morning cortisol response in patients with severe global amnesia; *Psychoneuroendocrinology* 2005 30, 101-105
- 20 Wessa et al. Altered cortisol awakening response in posttraumatic stress disorder; *Psychoneuroendocrinology* 2006 31, 209-215
- 21 Rohleder et al. Sex-specific adaptation of endocrine and inflammatory responses to repeated nauseogenic body rotation; *Psychoneuroendocrinology* 2006 31, 226-236
- 22 Rohleder et al. Stress on the dance floor: The cortisol stress response to social-evaluative threat in competitive ballroom dancers; *PSPB Vol. 33 No. 1, January* 2007 69-84

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LIABILITY: Complaints will be accepted in each mode –written or vocal. Preferred is that the complaint is accompanied with the test performance and results. Any modification of the test procedure or exchange or mixing of components of different lots could negatively affect the results. These cases invalidate any claim for replacement. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the kit during transportation is not subject to the liability of the manufacturer